# LUMINA™ IV dressing

en

STERILEEO

## Description:

LUMINA IV dressing consists of a thin film backing with a non-latex adhesive. The dressing is breathable, allowing good oxygen and moisture vapour exchange. It is waterproof and impermeable to liquids and bacteria. An intact dressing protects the site from outside contamination. The dressing has a slit adapted to be positioned around a catheter. The dressing should preferably be applied and removed by a healthcare professional. Before removal the adhesion is reduced by illuminating the dressing with the UV-A lamp or by exposing the dressing for sunlight.

## Indications:

LUMINA IV dressing is primarily intended to cover and protect catheter sites and wounds, and to secure devices to the skin.

## Contraindications:

Sensitivity towards UV-A radiation.  
Allergy / sensitivity towards acrylic adhesives.

## Warnings:

Do not use the dressing as a replacement for sutures and other primary wound closure methods.  
Do not look directly into the UV-A lamp when it is turned on.  
Do not expose the dressing to light during use, not even indoor light.  
If the sterile packaging is damaged, the dressing must not be used but discharged.

## Precautions:

1. Stop any bleeding at the site before applying the dressing.
2. Do not stretch the dressing during application as tension can cause skin trauma.
3. Make sure the skin is clean, free from soap residue and lotion, and allowed to dry thoroughly before applying the dressing to prevent skin irritation and to ensure good adhesion.
4. LUMINA IV dressings should not be re-sterilized.
5. Single use only.

## Instructions for use:

### Dressing selection:

Choose a dressing size large enough to provide at least two centimetres' margin which adheres to dry, healthy skin around the wound site.

### Site preparation:

1. Prepare the site according to institution protocol.
2. Clipping of hair at the site may improve dressing adhesion. Shaving may reduce pain on dressing removal.
3. Allow preps and protectants to dry completely before applying the dressing to prevent skin irritation and to ensure good adhesion.

### Application:

1. Open the package and take out the sterile dressing.
2. Peel off the middle part of the backing paper including securing strips and put aside.
3. Remove the release paper on the adhesive side.
4. Centre the dressing over the application site.
5. Press the dressing into place, and put the two flaps around the IV catheter.
6. Slowly remove the remaining part of the backing paper. Peel with a low angle to avoid lifting of the dressing.
7. Smooth the dressing from the centre toward the edges, using firm pressure to enhance adhesion.

8. Cover the dressing from direct or indirect sunlight with clothes or bed clothes, to prevent unintentional reduction of adhesion. (see pictures in Figure 1)

Figure 1

### Site care:

1. The site should be observed for signs of infection or other complications. If infection is suspected, remove the dressing, inspect the site directly, and determine appropriate medical intervention. Infection may be signalled by fever, pain, redness, swelling, or unusual odour or discharge.
2. Change the dressing according to institution protocol, or when the barrier properties have been compromised. For example, if the dressing is loose or punctured by sutures or other devices, it should be changed.

### Removal:

1. Take the UV-A lamp and direct it at the IV dressing at a distance of 3-5 centimetres (see pictures in Figure 2 below).
2. Turn on the lamp and illuminate an area of the IV dressing for 3-5 seconds.
3. Move the UV-A lamp to illuminate the next area of the IV dressing for 3-5 seconds.
4. Repeat step 3 until the entire area of the IV dressing has been illuminated (see picture 2).
5. Turn off the UV-A lamp.
6. Gently grasp an edge and slowly peel the IV dressing from the skin in the direction of hair growth (see picture 3). If adhesion level is still high, repeat step 1-6 with 5 seconds illumination per area.


Figure 2

### Shelf life and storage information:

For best results, store in a cool, dark and dry place. For shelf life, refer to the expiration date stated on each package.

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**Lumina Adhesives AB**